

## Test tubes with physiological saline solution



The test tubes with physiological saline solution NaCl 0.9% serve as transport vessels for swab samples.

The “ready to use” tubes are available in various sizes and designs and are supplied sterile. On request, the stoppers are also available in other colours as well as with individual (barcode) label.

Item no.	Description	Packing unit
042500	Test tubes 4 ml, 1 ml NaCl 0.9%	100 pieces/bag, 2,000 pieces/box
0961 1253	Test tubes 13 ml, 3 ml NaCl 0.9%	100 pieces/bag, 1,000 pieces/box
0961 3251	Test tubes 8 ml, 2 ml NaCl 0.9%	100 pieces/bag, 1,000 pieces/box
0961 3363	Test tubes 5 ml, 1,5 ml NaCl 0.9%	100 pieces/bag, 2,000 pieces/box
0961 325X	Test tubes 8 ml, 2 ml NaCl 0.9%, with individual label	On request

Upon request, we can also supply you with compatible swabs, as well as tubes whose screw caps are designed to hold the swab - which has been broken off and inserted into the tube. During automated processing, the swab is removed along with the tube as soon as it is opened, making the sample material easily accessible for processing. For sample storage, we offer suitable screw caps for resealing.



Item no.	Description	Packing unit
056504	Nylon flock fiber swab, nose breakpoint 85 mm, tip-Ø 2.5 mm	On request
056509	Nylon flock fiber swab, throat breakpoint 80 mm, tip-Ø 5 mm	On request